

CORRESPONDENCE

using halothane *via* the tracheostomy. Bronchoscopy revealed marked tracheomalacia and tracheal collapse during spontaneous ventilation. The tracheostomy was, therefore, left in place, and the child was subsequently discharged home without further adverse events.

To my knowledge, there are no reports of respiratory compromise secondary to oral midazolam administration. In addition, during the past 5 yr, I have used oral midazolam premedication in more than 100 children with sleep apnea presenting for tonsillectomy, and I have not seen worsening of upper airway obstruction. Respiratory compromise has been reported in healthy children after nasally administered midazolam.^{1,2} Peak plasma concentrations after intranasal midazolam are generally higher than those after oral administration.¹ This complication illustrates that, although oral midazolam premedication is generally extremely safe, selected patients with airway abnormalities might be at greater than usual risk for respiratory compromise.

Ronald S. Litman, D.O.

Assistant Professor of Anesthesiology and Pediatrics

Anesthesiology
1996;85:1218
© 1996 American Society of Anesthesiologists, Inc.
Lippincott-Raven Publishers

ParaGraph Muscle Stimulator: New Approach to Placement

To the Editor:—The ParaGraph Muscle Stimulator (Vital Signs, Totowa, NJ) is a new type of neuromuscular stimulator that uses a piezoelectric motion sensor. The sensor, when correctly applied, allows for the electronic interpretation and a computer screen display of the movement of a stimulated muscle. This permits monitoring of muscle relaxation without direct contact with the patient's extremity, a distinct advantage because an extremity is frequently not available during surgery.

The literature provided with the ParaGraph Muscle Stimulator* suggests placing the electric stimulator pad over the ulnar nerve at the wrist and the motion sensor over the bulk of the thenar eminence. We have found more satisfactory results by placing the motion sensor over the intersection of the hypothenar eminence and digit minimi (fifth digit). Stimulation of the ulnar nerve causes more motion of

* ParaGraph Operator's Manual. Vital Signs, Inc. 1993.

Anesthesiology
1996;85:1218-9
© 1996 American Society of Anesthesiologists, Inc.
Lippincott-Raven Publishers

In Reply:—As inventors and developers of the ParaGraph muscle relaxation monitor, we support the findings of Roberts and Dorsch regarding the placement of the piezo-electric motion sensor of the ParaGraph monitor.

Historically, the piezo-electric motion sensor was placed on the thumb to record movement of the adductor pollicis muscle to facilitate clinical comparison with an isometric mechanomyogram. The

Chief, Division of Pediatric Anesthesia
University of Rochester School of Medicine and Dentistry
Strong Memorial Hospital
601 Elmwood Avenue
Rochester, New York 14642
Electronic mail: Rlitman@ccmail.anes.rochester.edu

References

1. Malinovsky J, Populaire C, Cozian A, Lepage JY, Lejus C, Pinaud M: Premedication with midazolam in children. Effect of intranasal, rectal and oral routes on plasma midazolam concentrations. *Anaesthesia* 1995; 50:351-4.
2. Karl HW, Keifer AT, Rosenberger JL, Larach MG, Ruffe JM: Comparison of the safety and efficacy of intranasal midazolam or sufentanil for preinduction of anesthesia in pediatric patients. *ANESTHESIOLOGY* 1992; 76:209-15

(Accepted for publication August 20, 1996.)

the fifth digit than the thumb. Another approach is to place the stimulating pad more lateral, so the median nerve is stimulated. This will cause the thumb to move more.

It is important to visually assess the response to stimulation before administering muscle relaxants, to determine where the motion sensor and stimulating pad should be placed.

Christopher Roberts, M.D.

Resident in Anesthesia

Jerry A. Dorsch, M.D.

Consultant in Anesthesia

Department of Anesthesia and Critical Care
Mayo Clinic

4500 San Pablo Road

Jacksonville, Florida 32224

(Accepted for publication August 21, 1996.)

mechanomyogram is the clinical standard for measurement of muscle relaxation. It requires the hand be fixed to an arm board and the thumb to be placed in a ring so that force generated by the thumb in response to stimulation can be quantified. These research studies required placing the ParaGraph motion sensor over the thumb, to compare its response with the mechanomyogram. However, during routine clinical use, the motion sensor should be placed over the

COR
muscl
by Ro
We
Robe
1) Pla
tio
2) Ele
tw
3) O
in
4) Pl
wl
We
ment
given
tor. I
thema
Dors
Anest
1996;
© 199
Lippin
To
comp
froze
of er
sary,
that
froze
bloo
alwa
same
plas
regis
bloo
Anest
1996;
© 199
Lippin
In
freez
and t
able
tion
is the
Anest