

Date :

Title : INACTIVATION OF METHIONINE SYNTHETASE (MS) BY NITROUS OXIDE

Authors : D. D. Koblin, Ph.D., J. E. Watson, J. E. Deady, M.S., E. L. R. Stokstad, Ph.D., and E. I. Eger, II, M.D.

Affiliation: Department of Anesthesia, University of California, San Francisco, California 94143, and Department of Nutritional Sciences, University of California, Berkeley, California 94720

Introduction. Exposure of rats to 50% N₂O for 6 hr completely inactivates liver methionine synthetase (MS),¹ a vitamin B₁₂-dependent enzyme that synthesizes the amino acid methionine from homocysteine. The inactivation of this enzyme may play a role in the development of polyneuropathy and anemia seen in subjects chronically exposed to subanesthetic² and trace levels³ of nitrous oxide. Because the effects of this process may be important to patients and operating room personnel, we studied the time course of inactivation, its dependence on anesthetic concentration, the time course of recovery of activity, and the possibility that other anesthetics might also produce inactivation.

Methods. Groups of eight adult male CD-1 mice were placed in a 20-l hyperbaric chamber and exposed to a range of N₂O concentrations (0 to 0.8 atm N₂O) for varying periods of time (0 to 4 hr). They were also exposed to xenon (0.8 atm), nitrogen (0.8 atm), halothane (0.6% atm), isoflurane (0.6% atm), and enflurane (1.0% atm), all in the presence of 1 atm O₂. Anesthetic concentrations were determined by gas chromatography. Following anesthetic exposure, livers were removed and homogenized in 0.01 M phosphate buffer (pH 7.3). The homogenate was centrifuged at 20,000 g for 80 min. Methionine synthetase activity was measured by incubating an aliquot of the supernatant with [¹⁴CH₃]-tetrahydrofolate and homocysteine at 37 C for 1 hr, and then measuring the amount of radioactive methionine produced.⁴ Methionine synthetase activity is expressed as the nanomoles of methionine produced per hour per gram of liver.

Results. The MS activity in livers of mice exposed to 0.8 atm N₂O decreased progressively as exposure times increased (table 1). Periods of exposure as short as 15 to 30 min produced significant inactivation. Exposure to xenon did not affect MS activity (table 1). Methionine synthetase activity was progressively inhibited as N₂O concentrations increased (table 2). For 4-hr exposure periods, 50% inactivation was achieved by approximately 0.1 atm N₂O (table 2). Methionine synthetase activity returned to 70% of its control value one day after exposure to 0.8 atm N₂O for 4 hr (table 3). Recovery of enzyme activity was nearly complete two days after exposure to N₂O (table 3). Exposure of mice to anesthetizing levels of halothane, enflurane, and isoflurane did not affect MS activity significantly. We conclude that brief exposures to nitrous oxide inactivate an essential metabolic enzyme, an effect that appears to be confined to anesthetic or subanesthetic, but not trace, levels. Thus, our results may have greater implications for patients than for operating room personnel. Studying chronic exposure to trace levels may alter this impression of innocuousness to operating room personnel.

References

1. Deacon R, Lumb M, Perry J, et al: Selective inactivation of vitamin B₁₂ in rats by nitrous oxide. *Lancet* II:1023-1024, 1978
2. Amess JAL, Burman JF, Rees GM, et al: Megaloblastic hemopoiesis in patients receiving nitrous oxide. *Lancet* II:339-342, 1978
3. Cohen EN, Brown BW, Wu M, et al: Anesthetic health hazards in the dental operatory. *Anesthesiology* 51(3S):254, 1979 (abstract)
4. Sauer HJ, Jaenicke L: Einfacher Test zur Messung der Methionin-Synthetase-(MS-) Aktivität und seine Anwendungsmöglichkeiten in der Klinik. *Klin Wochenschr* 50:986-990, 1972

Table 1. Time Course of N₂O Inactivation

Group	MS Activity (mean ± SE) (nmol methionine/g/hr)
Shelf controls	296 ± 32.1
(15 min) 0.8 atm N ₂ O	197 ± 11.7
(30 min) 0.8 atm N ₂ O	54 ± 8.3
(60 min) 0.8 atm N ₂ O	56 ± 4.6
(120 min) 0.8 atm N ₂ O	41 ± 19.2
(240 min) 0.8 atm N ₂ O	16 ± 5.5
(240 min) 0.8 atm N ₂	166 ± 11.5
(240 min) 0.8 atm Xe	302 ± 30.7

Table 2. Concentration Dependence of N₂O Inactivation

Group	MS Activity (mean ± SE) (nmol methionine/g/hr)
Shelf controls	297 ± 68.3
0.003 atm N ₂ O (4 hr)	304 ± 18.9
0.010 atm N ₂ O (4 hr)	277 ± 5.9
0.049 atm N ₂ O (4 hr)	242 ± 10.7
0.20 atm N ₂ O (4 hr)	101 ± 3.1
0.79 atm N ₂ O (4 hr)	58 ± 3.0

Table 3. Recovery of MS Activity after a 4-Hr Exposure to 0.8 atm N₂O

Group	MS Activity (mean ± SE) (nmol methionine/g/hr)
Shelf controls	413 ± 32.0
0 day recovery	49 ± 9.0
1 day recovery	279 ± 14.2
2 day recovery	333 ± 25.0
4 day recovery	367 ± 40.1
8 day recovery	323 ± 14.8