

Fig. 1 Experimental Setup for PIV/LIF simultaneous measurement




Fig. 2 Different three mixing periods (A) : Continuous infusing period, (B) : Rapid turbulent mixing period, (C) : Slow molecular mixing period




Fig. 3 The temporal and spatial mixing characteristics in a stirred mixing tank

VISUALIZATION OF TURBULENT MASS TRANSFER IN A STIRRED MIXER

Eun Ho Jeong, Sang Youl Yoon and Kyung Chun Kim

Microfluidics and Bio-Fluids Laboratory, Pusan National Univ., Korea

Simultaneous measurements of turbulent velocity and concentration fields in a stirred mixer tank have been carried out by using a combined PIV/LIF technique. Instantaneous velocity fields are measured by a 1K×1K CCD camera while the concentration fields are determined by measuring the fluorescence intensity of Rhodamine B tracer excited by the second pulse of Nd:Yag laser light (Fig. 1). It is shown that the general features of the mixing

patterns are highly dependent on the local flow characteristics during the rapid decay of the mean concentration. However, the small scale mixing is found to be independent of the local turbulent velocity fluctuations. The correlation fields between concentration fluctuation and velocity fluctuation represents the active mixing region.